

# adapt

## Designed for adjustability

- Sustainable tilt in space mobility base, to accommodate multiple conditions and changing postural needs.
- Proven industry-wide seat system compatibility.
- Durability, longevity, reusability.



Discover More



## adapt

## Technical Data Sheet







Chair Dimensions and Sizing Guide

SPECIFICATION	XS	SMALL	Medium	Large	
Max User Mass	110kg	110kg	110kg	110kg	
Seat Cushion Depth Range	30-35cm	35-40cm	40-45cm	45-50cm	
Seat Width	30cm	35cm	40cm	45cm	
Back Rest Height	40cm	45cm	51cm	57cm	
Total Backrest Height Growth	5cm	5cm	5cm	5cm	
Backrest Width	25cm	30cm	31.5cm	35cm	
Leg Rest Angle Range	80°-170°	80°-170°	80°-170°	80°-170°	
Leg Rest Length Range	21-48cm	21-48cm	21-48cm	21-48cm	
Back Rest Recline Angle Range	90-180°	90-180°	90-180°	90-180°	
Tilt in Space	+3° to -45°	+3° to -45°	+3° to -45°	-45° +3° to -45°	



36-44cm

21-43cm

20-30cm

15-26cm

54-66cm

21-48cm

35-45cm

30-41cm

(14-17 in)

(8-17 in)

(8-12 in)

(6-10 in)


(21-26 in)

(8-19 in)

(14-18 in)

(12-16 in)





S			
42-52cm	(17-20in)		4
21-43cm	(8-17 in)		2
25-35cm	(10-14 in)		3
20-31cm	(8 -12 in)		2



## Enhancing access and independence for all

The foundation of manual wheelbase mobility: Chunc Adapt's modular design offers extensive customization options for clients who cannot maintain a symmetrical and unsupported sitting posture.

Adjustable components maintain the integrity of postural support while offering width and additional depth adjustments, allowing the system to grow with the existing user or be modified by a trained technician for a new client—delivering a lower lifetime cost.

- Indoor and outdoor use.



- Modular design, available in 4 sizes (XS, S, M, L) with growth adjustability.
- 45° Tilt Adjustment: support various seating positions, including hoisted transfers.
- Multi option mounting kit: fixed seat or rail mounting (15" & 18").
- Back rest inclination and width adjustment (up to 14").
- Options:
- Push bar, pram handle.
- Foot operated brake levers, drum brake assembly.
- Leg supports swing away & angle adjustable, individual, foot box.
- Compact & manoeuvrable, 1.2m turning circle.
- Max load capacity 110kg/24 stone.

- Optimal posture alignment, risk reduction, functional improvement and aesthetically pleasing design.
- 45° Tilt Adjustment; enhances comfort and accessibility, for individuals who may lack strength to change position.
- 45° change reduces pressure by up to 30%. Tilting at varying degrees throughout the day also provides relief.
- Back recline option up to 180°, to assist with personal care.
- Grows with your client and adapts to the evolving requirements of progressive conditions.
- Power assist option confidence and freedom to do more.





## #thechuncway

Custom-made indoor and outdoor wheelchairs with built-in postural support for children and adults with moderate to complex needs.

Designed, manufactured and hand-assembled in Herefordshire, England (in 10 working days\*).

Supporting your everyday needs – and passionate about it too.

\*97% all standard chairs made in 2023



**Discover More** 

Web: chunc.com | Tel: 01432 377 512 | Email: info@chunc.co.uk